



Special Female Foley Catheter.
The design follows the anatomy and physiology of female urethra.

Silicone Foley Catheter – Female Type



Bioptimal™

www.bioptimalg.com

PCT Patent:

PCT/CN2018/073043

PCT/CN2018/073316

★ Characteristics

- 1 Circumferential positioning
- 2 Stably low Insertion
- 3 Complete Isolation
- 4 Uniform Intubation
- 5 Precise Guidance

★ Clinical expectation

To reduce the risk of microbe taken into urethra and bladder while inserting a catheter. To reduce injury to urethral mucosa caused by catheterization.

★ Preferred departments

Applicable to female patients in all clinical departments.

Model No.	Description	Size (FR/Ch)	Balloon Size (mL)	Length(mm)
FC214W I	2-way	14	5-15	410
FC214W II	2-way	14	20-30	410
FC216W I	2-way	16	5-15	410
FC216W II	2-way	16	20-30	410
FC218W I	2-way	18	5-15	410
FC218W II	2-way	18	20-30	410
FC220W I	2-way	20	5-15	410
FC220W II	2-way	20	20-30	410
FC314W I	3-way	14	5-15	410
FC314W II	3-way	14	20-30	410
FC316W I	3-way	16	5-15	410
FC316W II	3-way	16	20-30	410
FC318W I	3-way	18	5-15	410
FC318W II	3-way	18	20-30	410
FC320W I	3-way	20	5-15	410
FC320W II	3-way	20	20-30	410

BIOPTIMAL INTERNATIONAL PTE. LTD.

36 Jalan Tukang
Singapore 619266

Tel: +65 6213 5777

Fax: +65 6213 5737

Email: sales@bioptimalg.com

