

Vascular Introducer Kit

prUNUS



Provider of Critical Care Product

PRODUCT DESCRIPTION

Vascular Introducer Kit is a basic customized kit that comprises of an introducer with or without a vessel dilator, a guidewire, together with some other components mainly to assist in the catheterization procedures.

PRINCIPLE OF OPERATION

Vascular Introducer Kit is intended to assist in the introduction of temporary pacing leads, balloon catheters, closed end catheter or angiographic catheters into the vasculature.

SPECIFICATIONS

Model No.	Size (Fr.)	Length (cm)	Guidewire
CPIK6F	6F	10cm	0.035" (0.89mm) × 60cm
CPIK7F	7F	10cm	0.035" (0.89mm) × 60cm
CPIK8F	8F	10cm	0.035" (0.89mm) × 60cm
CPIK8.5F	8.5F	10cm	0.035" (0.89mm) × 60cm

Optional Accessories

18G x 2 1/2" Needle

Syringe Luer Slip 5mL/10mL

OTN Catheter

Guide Wire (0.035"/0.89mm- J60cm)

Scalpel # 11-Short

Gauze Swab 7.5x7.5x8PLY

3-0 Sutures

Paper Towel 17" x 22"

Needle 25G or 22G

Drape Minor Proc. Fen 22x22"

Syringe 1.5cc

Note: Customization is available, please contact your local sales representative for more information.

Bioptimal International Pte. Ltd.

36 Jalan Tukang, Singapore 619266
Tel: +65 6213 5777 Email: Sales@bioptimalg.com
Web: <https://www.bioptimalg.com>
Version: BRH-00-01 Rev.02